

**APRIL 11, 1997 VA HANDBOOK 5383.2**  
**APPENDIX A**  
**COLLECTION SITE SUPPLIES**

**1. Blueing Material.** In order to avoid staining the toilet bowl or urinal, it is recommended that the blueing material be a chemical made for use in such fixtures.

Suggested product: Blue VANISH (Liquid)

Vendor: Drackett Company

5020 Spring Grove Avenue

Cincinnati, OH 45232

Telephone: (513) 632-1500

**2. Collection Kit**

- a. Bottle A: Pre-sealed 100 mL plastic urine specimen bottle with affixed Temperature Monitoring Thermometer.
- b. Bottle B: Pre-sealed 100 mL plastic urine specimen bottle.
- c. Gray foam insert to secure urine specimen bottles.
- d. Specimen Bag, 8" X 8" leak proof ziploc bag containing one liquid absorbing sheet.
- e. Mailing box with kit integrity seal.
- f. Urine Specimen Collection Cup (Wide Mouth Collection Cup - for females), 105 mL Urine Specimen Collection Cup, pre-sealed in integrity pouches.
- g. Federal Drug Administration (FDA) Insert, printed on Blue stock with Black Ink. Information is included as a requirement of FDA.
- h. Kit will have a pre-printed mailing label affixed to the box and addressed to the testing laboratory.
- i. Two serrated red strips to cover bottles.

**3. Oral Thermometer.** Preferably digital with disposable tips.

**4. Chain of Custody Forms**

- a. OMB Form No. 9999-0023 (form approved for Government-wide use).
- b. Replenishment available from the VA Service and Distribution Center (Stock Number F06088).

**5. Miscellaneous Items**

- a. Disposable drinking glasses.
- b. Writing implements, e.g., pens, pencils, etc.
- c. Folder to retain laboratory copies of chain of custody forms.